


駆動部品名



※は受注発注品です。

駆動部品名	バキューミキサー					バキュー スパット
	#6400 984-226	#4450 984-227	#6500 984-228	#6600 984-230	#7600(金属製) 984-231	#4120 984-221
	200mL	300mL	500mL	875mL	1200mL	200mL
						
	φ85 (mm)	φ85 (mm)	φ105 (mm)	φ105 (mm)	φ120 (mm)	φ85 (mm)
①ドライブナット	#6375(低速用) 984-248				※#6865	#4109(高速用) 984-246
②リド(フタ)	#6401 984-267		#6501 984-268		※#7601 984-272	#4001 984-259
③ボウル	#4002 984-261	#4452 984-262	#6502 984-263	#6602 984-264	※#7602 984-266	#4002 984-261
④パドル(羽根)	#6406(石膏用) 984-274	#4455(石膏用) 984-276	#6506(石膏用) 984-277 #6508 (アルギン酸印象材用) 984-279	#6606 (アルギン酸印象材用) 984-280	※#7606 984-282	#4007 984-273

製品名：バキューベスター コンビネーション ユニット 分類：一般医療機器 医療機器届出番号：13B1X00004000033 一般の名称：歯科技工用真空攪拌器